

## 兔抗 SHPRH 多克隆抗体

中文名称：兔抗 SHPRH 多克隆抗体

英文名称： Anti-SHPRH rabbit polyclonal antibody

别名： bA545I5.2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SHPRH

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	SHPRH is a ubiquitously expressed protein that contains motifs characteristic of several DNA repair proteins, transcription factors, and helicases. SHPRH is a functional homolog of <i>S. cerevisiae</i> RAD5 (Unk et al., 2006 [PubMed 17108083]).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SHPRH
<b>Immunogen:</b>	Fusion protein of human SHPRH
<b>Full name:</b>	SNF2 histone linker PHD RING helicase, E3 ubiquitin protein ligase
<b>Synonyms :</b>	bA545I5.2

<b>SwissProt:</b>	Q149N8
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	30-150

