

兔抗 SHROOM2 多克隆抗体

中文名称：兔抗 SHROOM2 多克隆抗体

英文名称： Anti-SHROOM2 rabbit polyclonal antibody

别名： APXL; HSAPXL

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： SHROOM2

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The protein encoded by this gene shares significant similarities with the apical protein from *Xenopus laevis* which is implicated in amiloride-sensitive sodium channel activity. This gene is a strong candidate gene for oc

	ular albinism type 1 syndrome.
Applications:	ELISA, IHC
Name of antibody:	SHROOM2
Immunogen:	Synthetic peptide of human SHROOM2
Full name:	shroom family member 2
Synonyms :	APXL; HSAPXL
SwissProt:	Q13796
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer and Human cervical cancer
IHC Recommend dilution:	50-200

