

## 兔抗 SIAH2 多克隆抗体

中文名称：兔抗 SIAH2 多克隆抗体

英文名称： Anti-SIAH2 rabbit polyclonal antibody

别名： hSiah2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SIAH2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

|                          |  |
|--------------------------|--|
| <b>Background:</b>       | This gene encodes a protein that is a member of the seven in absentia homolog (SIAH) family. The protein is an E3 ligase and is involved in ubiquitination and proteasome-mediated degradation of specific proteins. The activity of this ubiquitin ligase has been implicated in regulating cellular response to hypoxia. |
| <b>Applications:</b>     | ELISA, IHC   |
| <b>Name of antibody:</b> | SIAH2  |
| <b>Immunogen:</b>        | Synthetic peptide of human SIAH2   |
| <b>Full name:</b>        | siah E3 ubiquitin protein ligase 2   |
| <b>Synonyms :</b>        | hSiah2   |

|                                    |                                     |
|------------------------------------|-------------------------------------|
| <b>SwissProt:</b>                  | O43255                              |
| <b>ELISA Recommended dilution:</b> | 1000-2000                           |
| <b>IHC positive control:</b>       | Human brain and Human breast cancer |
| <b>IHC Recommend dilution:</b>     | 10-50                               |

