

## 兔抗 SIGLEC6 多克隆抗体

中文名称：兔抗 SIGLEC6 多克隆抗体

英文名称：Anti-SIGLEC6 rabbit polyclonal antibody

别名：sialic acid binding Ig like lectin 6; CD327; CD33L; OBBP1; CD33L1; CD33L2; CDW327

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Rabbit

抗原：SIGLEC6

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the SIGLEC (sialic acid binding immunoglobulin-like lectin) family of proteins. The encoded transmembrane receptor binds sialyl-TN glycans and leptin. Placental expression of the encoded protein is upregulated in preeclampsia.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SIGLEC6
<b>Immunogen:</b>	Fusion protein of human SIGLEC6
<b>Full name:</b>	sialic acid binding Ig like lectin 6
<b>Synonyms:</b>	CD327; CD33L; OBBP1; CD33L1; CD33L2; CDW327

<b>SwissProt:</b>	O43699
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human breast cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	40-200
<b>WB Predicted band size:</b>	50 kDa
<b>WB Positive control:</b>	Human right lower lung tissue and Human plasma s olution lysates
<b>WB Recommended dilution:</b>	500-2000



