

## 兔抗 SIGLEC9 多克隆抗体

中文名称：兔抗 SIGLEC9 多克隆抗体

英文名称： Anti-SIGLEC9 rabbit polyclonal antibody

别名： CD329; CDw329; FOAP-9; siglec-9; OBBP-LIKE

相关类别： 一抗

抗原： SIGLEC9

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Sialic acid-binding Ig-like lectin 9 is a protein that in humans is encoded by the SIGLEC9 gene. Putative adhesion molecule that mediates sialic-acid dependent binding to cells. Preferentially binds to alpha-2,3- or alpha-2,6-linked sialic acid. The sialic acid recognition site may be masked by cis interactions with sialic acids on the same cell surface. Expressed by peripheral blood leukocytes (neutrophils and monocytes but not eosinophils).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SIGLEC9
<b>Immunogen:</b>	Fusion protein of human SIGLEC9

<b>Full name:</b>	sialic acid binding Ig-like lectin 9
<b>Synonyms :</b>	CD329; CDw329; FOAP-9; siglec-9; OBBP-LIKE
<b>SwissProt:</b>	Q9Y336
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and human brain
<b>IHC Recommend dilution:</b>	50-200

