

## 兔抗 SIGLEC10 多克隆抗体

中文名称：兔抗 SIGLEC10 多克隆抗体

英文名称： Anti-SIGLEC10 rabbit polyclonal antibody

别名： SLG2; PRO940; SIGLEC-10

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SIGLEC10

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	SIGLECs are members of the immunoglobulin super family that are expressed on the cell surface. Most SIGLECs have 1 or more cytoplasmic immune receptor tyrosine-based inhibitory motifs, or ITIMs. SIGLECs are typically expressed on cells of the innate immune system, with the exception of the B-cell expressed SIGLEC6.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SIGLEC10
<b>Immunogen:</b>	Synthetic peptide of human SIGLEC10
<b>Full name:</b>	sialic acid binding Ig-like lectin 10

<b>Synonyms :</b>	SLG2; PRO940; SIGLEC-10
<b>SwissProt:</b>	Q96LC7
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	25-100

