

兔抗 XPO1 多克隆抗体

中文名称：兔抗 XPO1 多克隆抗体

英文名称： Anti-XPO1 rabbit polyclonal antibody

别名： exportin 1; emb; CRM1; exp1

相关类别： 一抗

抗原： XPO1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This cell-cycle-regulated gene encodes a protein that mediates leucine-rich nuclear export signal (NES)-dependent protein transport. The protein specifically inhibits the nuclear export of Rev and U snRNAs. It is involved in the control of several cellular processes by controlling the localization of cyclin B, MPAK, and MAPKAP kinase 2. This protein also regulates NFAT and AP-1.
Applications:	ELISA, WB, IHC
Name of antibody:	XPO1
Immunogen:	Synthetic peptide of human XPO1
Full name:	exportin 1

Synonyms:	emb; CRM1; exp1
SwissProt:	O14980
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer and human esophagus cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	123 kDa
WB Positive control:	Hela,K562 and HEPG2 cell lysates
WB Recommended dilution:	500-2000



