

兔抗 XPO6 多克隆抗体

中文名称：兔抗 XPO6 多克隆抗体

英文名称：Anti-XPO6 rabbit polyclonal antibody

别名：EXP6; RANBP20

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：XPO6

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the importin-beta family. Members of this family are regulated by the GTPase Ran to mediate transport of cargo across the nuclear envelope. This protein has been shown to mediate nuclear export of profilin-actin complexes. A pseudogene of this gene is located on the long arm of chromosome 14. Alternative splicing results in multiple transcript variants that encode different protein isoforms.
Applications:	ELISA, IHC
Name of antibody:	XPO6
Immunogen:	Fusion protein of human XPO6

Full name:	exportin 6
Synonyms:	EXP6; RANBP20
SwissProt:	Q96QU8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and human ovarian cancer
IHC Recommend dilution:	25-100



