

## 兔抗 PLA1A 多克隆抗体

中文名称: 兔抗 PLA1A 多克隆抗体

英文名称: Anti-PLA1A rabbit polyclonal antibody

别名: PSPLA1; PS-PLA1

储存: 冷冻 (-20℃) 避光

抗原: PLA1A

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a phospholipase that hydrolyzes fatty acids at the sn-1 position of phosphatidylserine and 1-acyl-2-lysophosphatidylserine. This secreted protein hydrolyzes phosphatidylserine in liposomes. Three transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	IHC
<b>Name of antibody:</b>	PLA1A
<b>Immunogen:</b>	Synthesized peptide derived from internal of human PLA1A.
<b>Full name:</b>	phospholipase A1 member A
<b>Synonyms :</b>	PSPLA1; PS-PLA1
<b>SwissProt:</b>	Q53H76

<b>IHC positive control:</b>	Human brain tissue
<b>IHC Recommend dilution:</b>	50-100

