

兔抗 PLA2G12A 多克隆抗体

中文名称： 兔抗 PLA2G12A 多克隆抗体

英文名称： Anti-PLA2G12A rabbit polyclonal antibody

别名： GXII; ROSSY; PLA2G12

相关类别： 一抗

抗原： PLA2G12A

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Secreted phospholipase A2 (sPLA2) enzymes liberate arachidonic acid from phospholipids for production of eicosanoids and exert a variety of physiologic and pathologic effects. Group XII sPLA2s, such as PLA2G12A, have relatively low specific activity and are structurally and functionally distinct from other sPLA2s (Gelb et al., 2000 [PubMed 11031251]).
Applications:	ELISA, IHC
Name of antibody:	PLA2G12A
Immunogen:	Fusion protein of human PLA2G12A
Full name:	phospholipase A2, group XIIA

Synonyms :	GXII; ROSSY; PLA2G12
SwissProt:	Q9BZM1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	20-100

