

兔抗 XRCC2 多克隆抗体

中文名称：兔抗 XRCC2 多克隆抗体

英文名称：Anti-XRCC2 rabbit polyclonal antibody

别名：DKFZp781P091

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：XRCC2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the RecA/Rad51-related protein family that participates in homologous recombination to maintain chromosome stability and repair DNA damage. This gene is involved in the repair of DNA double-strand breaks by homologous recombination and it functionally complements Chinese hamster irs 1, a repair-deficient mutant that exhibits hypersensitivity to a number of different DNA-damaging agents.
Applications:	WB
Name of antibody:	XRCC2
Immunogen:	Fusion protein of human XRCC2

Full name:	X-ray repair complementing defective repair in Chinese hamster cells 2
Synonyms :	DKFZp781P091
SwissProt:	O43543
WB Predicted band size:	32 kDa
WB Positive control:	HeLa, SH-SY5Y and SW480 cells
WB Recommended dilution:	500-2000

