

## 兔抗 XRCC3 多克隆抗体

中文名称：兔抗 XRCC3 多克隆抗体

英文名称： Anti-XRCC3 rabbit polyclonal antibody

别名： CMM6

抗原： XRCC3

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a member of the RecA/Rad51-related protein family that participates in homologous recombination to maintain chromosome stability and repair DNA damage. This gene functionally complements Chinese hamster irs1SF, a repair-deficient mutant that exhibits hypersensitivity to a number of different DNA-damaging agents and is chromosomally unstable. A rare microsatellite polymorphism in this gene is associated with cancer in patients of varying radio sensitivity.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	XRCC3
<b>Immunogen:</b>	Synthesized peptide derived from Internal of human

	XRCC3.
<b>Full name:</b>	X-ray repair complementing defective repair in Chinese hamster cells 3
<b>Synonyms :</b>	CMM6
<b>SwissProt:</b>	O43542
<b>IHC positive control:</b>	Human lung carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	38 kDa
<b>WB Positive control:</b>	HepG2 cells treated with Adriamycin (0.5uM, 5hours) lysate
<b>WB Recommended dilution:</b>	500-3000



