

兔抗 SIPA1 多克隆抗体

中文名称：兔抗 SIPA1 多克隆抗体

英文名称：Anti-SIPA1 rabbit polyclonal antibody

别名：SPA1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SIPA1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The product of this gene is a mitogen induced GTPase activating protein (GAP). It exhibits a specific GAP activity for Ras-related regulatory proteins Rap1 and Rap2, but not for Ran or other small GTPases. This protein may also hamper mitogen-induced cell cycle progression when abnormally or prematurely expressed. It is localized to the perinuclear region. Two alternatively spliced variants encoding the same isoform have been characterized to date.
Applications:	ELISA, IHC
Name of antibody:	SIPA1

Immunogen:	Synthetic peptide of human SIPA1
Full name:	signal-induced proliferation-associated 1
Synonyms :	SPA1
SwissProt:	Q96FS4
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer and Human cervical cancer
IHC Recommend dilution:	50-200

