

兔抗 XRCC5 多克隆抗体

中文名称：兔抗 XRCC5 多克隆抗体

英文名称： Anti-XRCC5 rabbit polyclonal antibody

别名： KU80; KUB2; Ku86; NFIV; KARP1; KARP-1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： XRCC5

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is the 80-kilodalton subunit of the Ku heterodimer protein which is also known as ATP-dependant DNA helicase II or DNA repair protein XRCC5. Ku is the DNA-binding component of the DNA-dependent protein kinase, and it functions together with the DNA ligase IV-XRCC4 complex in the repair of DNA double-strand break by non-homologous end joining and the completion of V(D)J recombination events.
Applications:	WB, IHC
Name of antibody:	XRCC5
Immunogen:	Synthesized peptide derived from C-terminal of human XRCC5.
Full name:	X-ray repair complementing defective repair in Chinese ha

	mster cells 5 (double-strand-break rejoining)
Synonyms :	KU80; KUB2; Ku86; NFIV; KARP1; KARP-1
SwissProt:	P13010
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	83 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000



