

## 兔抗 XRCC6 多克隆抗体

中文名称：兔抗 XRCC6 多克隆抗体

英文名称： Anti-XRCC6 rabbit polyclonal antibody

别名： ML8, KU70, TLA, CTC75, CTCBF, G22P1

抗原： XRCC6

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The p70/p80 autoantigen is a nuclear complex consisting of two subunits with molecular masses of approximately 70 and 80 kDa. The complex functions as a single-stranded DNA-dependent ATP-dependent helicase. The complex may be involved in the repair of nonhomologous DNA ends such as that required for double-strand break repair, transposition, and V(D)J recombination. High levels of autoantibodies to p70 and p80 have been found in some patients with systemic lupus erythematosus.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	XRCC6
<b>Immunogen:</b>	Fusion protein of human XRCC6

<b>Full name:</b>	X-ray repair complementing defective repair in Chinese hamster cells 6
<b>Synonyms :</b>	ML8, KU70, TLAA, CTC75, CTCBF, G22P1
<b>SwissProt:</b>	P12956
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human colon cancer and human ovarian cancer
<b>IHC Recommend dilution:</b>	10-50

