

## 兔抗 PLA2G4D 多克隆抗体

中文名称：兔抗 PLA2G4D 多克隆抗体

英文名称：Anti-PLA2G4D rabbit polyclonal antibody

别名：cPLA2delta

相关类别：一抗

抗原：PLA2G4D

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The phospholipase A2 enzyme family, including PLA2G4 D, catalyze the hydrolysis of glycerophospholipids at the sn-2 position and then liberate free fatty acids and lyso phospholipids (Chiba et al., 2004 [PubMed 14709560]).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PLA2G4D
<b>Immunogen:</b>	Fusion protein of human PLA2G4D
<b>Full name:</b>	phospholipase A2, group IVD (cytosolic)
<b>Synonyms :</b>	cPLA2delta
<b>SwissProt:</b>	Q86XP0
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human ovarian cancer
<b>IHC Recommend dilution:</b>	25-100

