

兔抗 SIRT4 多克隆抗体

中文名称：兔抗 SIRT4 多克隆抗体

英文名称： Anti-SIRT4 rabbit polyclonal antibody

别名： SIR2L4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SIRT4

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the sirtuin family of proteins, homologs to the yeast Sir2 protein. Members of the sirtuin family are characterized by a sirtuin core domain and grouped into four classes. The functions of human sirtuins have not yet been determined; however, yeast sirtuin proteins are known to regulate epigenetic gene silencing and suppress recombination of rDNA. Studies suggest that the human sirtuins may function as intracellular regulatory proteins with mono-ADP-ribosyltransferase activity.
Applications:	ELISA, IHC

Name of antibody:	SIRT4
Immunogen:	Synthetic peptide of human SIRT4
Full name:	sirtuin 4
Synonyms :	SIR2L4
SwissProt:	Q9Y6E7
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and Human lung cancer
IHC Recommend dilution:	25-100

