

## 兔抗 SIX5 多克隆抗体

中文名称：兔抗 SIX5 多克隆抗体

英文名称： Anti-SIX5 rabbit polyclonal antibody

别名： BOR2; DMAHP

抗原： SIX5

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a homeodomain-containing transcription factor that appears to function in the regulation of organogenesis. This gene is located downstream of the dystrophia myotonica-protein kinase gene. Mutations in this gene are a cause of branchiootorenal syndrome type 2.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	SIX5
<b>Immunogen:</b>	Synthesized peptide derived from internal of human SIX5.
<b>Full name:</b>	SIX homeobox 5

<b>Synonyms :</b>	BOR2; DMAHP
<b>SwissProt:</b>	Q8N196
<b>WB Predicted band size:</b>	75 kDa
<b>WB Positive control:</b>	K562 cells lysate
<b>WB Recommended dilution:</b>	500-3000

