

兔抗 YAP1 多克隆抗体

- 中文名称：兔抗 YAP1 多克隆抗体
- 英文名称：Anti-YAP1 rabbit polyclonal antibody
- 别名：Yes associated protein 1 YAP; YKI; COB1; YAP2; YAP65
- 相关类别：一抗
- 抗原：YAP1
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a downstream nuclear effector of the Hippo signaling pathway which is involved in development, growth, repair, and homeostasis. This gene is known to play a role in the development and progression of multiple cancers as a transcriptional regulator of this signaling pathway and may function as a potential target for cancer treatment. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	YAP1
Immunogen:	Synthetic peptide of human YAP1

Full name:	Yes associated protein 1
Synonyms:	YAP; YKI; COB1; YAP2; YAP65
SwissProt:	P46937
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human gastric cancer; Human ovarian cancer
IHC Recommend dilution:	20-100



