

兔抗 PLAAT4 多克隆抗体

中文名称：兔抗 PLAAT4 多克隆抗体

英文名称： Anti-PLAAT4 rabbit polyclonal antibody

别名： phospholipase A and acyltransferase 4; RIG1; TIG3;

抗原： PLAAT4

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Retinoids exert biologic effects such as potent growth inhibitory and cell differentiation activities and are used in the treatment of hyperproliferative dermatological diseases. These effects are mediated by specific nuclear receptor proteins that are members of the steroid and thyroid hormone receptor superfamily of transcriptional regulators. RARRES1, RARRES2, and RARRES3 are genes whose expression is upregulated by the synthetic retinoid tazarotene. RARRES3 is thought to act as a tumor suppressor or growth regulator. [provided by RefSeq, Jul 2008]
Applications:	ELISA, IHC

Name of antibody:	PLAAT4
Immunogen:	Fusion protein of human PLAAT4
Full name:	phospholipase A and acyltransferase 4
Synonyms:	RIG1; TIG3; HRSL4; HRASLS4; PLAAT-4; RARRES3; PLA1/2-3
SwissProt:	Q9UL19
IHC positive control:	Human colon cancer
IHC Recommend dilution:	25-100

