

兔抗 YARS2 多克隆抗体

中文名称：兔抗 YARS2 多克隆抗体

英文名称： Anti-YARS2 rabbit polyclonal antibody

别名： TYRRS; CGI-04; MLASA2; MT-TYRRS

相关类别： 一抗

抗原： YARS2

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a mitochondrial protein that catalyzes the attachment of tyrosine to tRNA(Tyr). Mutations in this gene are associated with myopathy with lactic acidosis and sideroblastic anemia type 2 (MLASA 2). Catalyzes the attachment of tyrosine to tRNA(Tyr) in a two-step reaction: tyrosine is first activated by ATP to form Tyr-AMP and then transferred to the acceptor end of tRNA(Tyr).
Applications:	ELISA, IHC
Name of antibody:	YARS2
Immunogen:	Fusion protein of human YARS2
Full name:	tyrosyl-tRNA synthetase 2, mitochondrial
Synonyms :	TYRRS; CGI-04; MLASA2; MT-TYRRS

SwissProt:	Q9Y2Z4
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human colon cancer and human lung cancer
IHC Recommend dilution:	50-200

