

## 兔抗 KCNS2 多克隆抗体

中文名称: 兔抗 KCNS2 多克隆抗体

英文名称: Anti-KCNS2 rabbit polyclonal antibody

别名: KV9.2

抗原: KCNS2

储存: 冷冻 (-20°C) 避光

宿主: Rabbit

反应种属: Human Mouse Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

|                          |   |
|--------------------------|---|
| <b>Background:</b>       | Potassium channel subunit. Modulates channel activity and reduces the ion flow By similarity. |
| <b>Applications:</b>     | WB  |
| <b>Name of antibody:</b> | KCNS2   |
| <b>Immunogen:</b>        | Synthesized peptide derived from internal of human KCNS2.                                     |
| <b>Full name:</b>        | potassium voltage-gated channel, delayed-rectifier, subfamily S, member 2                     |
| <b>Synonyms :</b>        | KV9.2   |
| <b>SwissProt:</b>        | Q9ULS6  |

|                                 |                    |
|---------------------------------|--------------------|
| <b>WB Predicted band size:</b>  | 54 kDa             |
| <b>WB Positive control:</b>     | HepG2 cells lysate |
| <b>WB Recommended dilution:</b> | 500-3000           |

