

兔抗 SIX6 多克隆抗体

中文名称: 兔抗 SIX6 多克隆抗体

英文名称: Anti-SIX6 rabbit polyclonal antibody

别名: SIX homeobox 6; Six9; ODRMD; OPTX2; MCOPCT2

抗原: SIX6

储存: 冷冻 (-20°C)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a homeobox protein that is similar to the Drosophila 'sine oculis' gene product. This gene is found in a cluster of related genes on chromosome 14 and is thought to be involved in eye development. Defects in this gene are a cause of isolated microphthalmia with cataract type 2 (MCOPCT2).
Applications:	ELISA, IHC
Name of antibody:	SIX6
Immunogen:	Synthetic peptide of human SIX6
Full name:	SIX homeobox 6

Synonyms:	Six9; ODRMD; OPTX2; MCOPCT2
SwissProt:	O95475
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	30-150

