

## 兔抗 SKI 多克隆抗体

中文名称：兔抗 SKI 多克隆抗体

英文名称：Anti-SKI rabbit polyclonal antibody

别名：C-SKI

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SKI

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	May play a role in terminal differentiation of skeletal muscle cells but not in the determination of cells to the myogenic lineage. Functions as a repressor of TGF-beta signaling.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SKI
<b>Immunogen:</b>	Synthetic peptide of humanSKI
<b>Full name:</b>	v-ski avian sarcoma viral oncogene homolog
<b>Synonyms :</b>	C-SKI
<b>SwissProt:</b>	P12755
<b>ELISA Recommended dilution:</b>	2000-10000

<b>IHC positive control:</b>	Human thyroid cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	100-300

