

兔抗 SLAMF1 多克隆抗体

中文名称：兔抗 SLAMF1 多克隆抗体

英文名称： Anti-SLAMF1 rabbit polyclonal antibody

别名： SLAM; CD150; CDw150

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLAMF1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	High-affinity self-ligand important in bidirectional T-cell to B-cell stimulation. SLAM-induced signal-transduction events in T-lymphocytes are different from those in B-cells. Two modes of SLAM signaling are likely to exist: one in which the inhibitor SH2D1A acts as a negative regulator and another in which protein-tyrosine phosphatase 2C (PTPN11)-dependent signal transduction operates.
Applications:	ELISA, WB
Name of antibody:	SLAMF1
Immunogen:	Fusion protein of human SLAMF1

Full name:	signaling lymphocytic activation molecule family member 1
Synonyms :	SLAM; CD150; CDw150
SwissProt:	Q13291
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	37 kDa
WB Positive control:	Jurkat cells
WB Recommended dilution:	200-1000

