

兔抗 PLBD1 多克隆抗体

中文名称：兔抗 PLBD1 多克隆抗体

英文名称：Anti-PLBD1 rabbit polyclonal antibody

别名：phospholipase B domain containing 1

相关类别：一抗

抗原：PLBD1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	In view of the small size of the putative binding pocket, it has been proposed that it may act as an amidase or a peptidase (By similarity). Exhibits a weak phospholipase activity, acting on various phospholipids, including phosphatidylcholine, phosphatidylinositol, phosphatidylethanolamine and lysophospholipids.
Applications:	ELISA, IHC
Name of antibody:	PLBD1
Immunogen:	Synthetic peptide of human PLBD1
Full name:	phospholipase B domain containing 1
SwissProt:	Q6P4A8

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	30-150

