

兔抗 YEATS2 多克隆抗体

中文名称： 兔抗 YEATS2 多克隆抗体

英文名称： Anti-YEATS2 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

抗 原： YEATS2

反应种属： Human

标 记 物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	YEATS2 is a scaffolding subunit of the ADA2A (TADA2A; MI M 602276)-containing (ATAC) histone acetyltransferase complex (Wang et al., 2008)
Applications:	WB, IHC
Name of antibody:	YEATS2
Immunogen:	Synthesized peptide derived from internal of human YEATS2.
Full name:	YEATS domain containing 2
Synonyms :	0
SwissProt:	Q9ULM3
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	151 kDa
WB Positive control:	293 cells lysate

WB Recommended dilution: 500-3000

