

兔抗 SLAMF8 多克隆抗体

中文名称：兔抗 SLAMF8 多克隆抗体

英文名称： Anti-SLAMF8 rabbit polyclonal antibody

别名： BLAME; CD353; SBBI42

相关类别： 一抗

抗原： SLAMF8

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the CD2 family of cell surface proteins involved in lymphocyte activation. These proteins are characterized by Ig domains. This protein is expressed in lymphoid tissues, and studies of a similar protein in mouse suggest that it may function during B cell lineage commitment. The gene is found in a region of chromosome 1 containing many CD2 genes.
Applications:	ELISA, IHC
Name of antibody:	SLAMF8
Immunogen:	Fusion protein of human SLAMF8

Full name:	SLAM family member 8p
Synonyms :	BLAME; CD353; SBBI42
SwissProt:	Q9P0V8
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human liver cancer and human thyroid cancer
IHC Recommend dilution:	50-200

