

兔抗 YES1 多克隆抗体

中文名称： 兔抗 YES1 多克隆抗体

英文名称： Anti-YES1 rabbit polyclonal antibody

别名： Yes; c-yes; HsT441; P61-YES

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： YES1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is the cellular homolog of the Yamaguchi sarcoma virus oncogene. The encoded protein has tyrosine kinase activity and belongs to the src family of proteins. This gene lies in close proximity to thymidylate synthase gene on chromosome 18, and a corresponding pseudogene has been found on chromosome 22.
Applications:	ELISA, WB
Name of antibody:	YES1
Immunogen:	Synthetic peptide of human YES1
Full name:	v-yes-1 Yamaguchi sarcoma viral oncogene homolog 1
Synonyms :	Yes; c-yes; HsT441; P61-YES

SwissProt:	P07947
ELISA Recommended dilution:	2000-10000
WB Predicted band size:	61 kDa
WB Positive control:	Lncap cells
WB Recommended dilution:	1000-5000

