

兔抗 SLC1A3 多克隆抗体

中文名称：兔抗 SLC1A3 多克隆抗体

英文名称： Anti-SLC1A3 rabbit polyclonal antibody

别名： EA6; EAAT1; GLAST; GLAST1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC1A3

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of a member of a high affinity glutamate transporter family. This gene functions in the termination of excitatory neurotransmission in central nervous system. Mutations are associated with episodic ataxia, Type 6. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB
Name of antibody:	SLC1A3
Immunogen:	Synthetic peptide of human SLC1A3
Full name:	solute carrier family 1 (glial high affinity glutamate transporter), member 3

Synonyms :	EA6; EAAT1; GLAST; GLAST1
SwissProt:	P43003
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	60 kDa
WB Positive control:	Human fetal brain tissue
WB Recommended dilution:	500-2000

