

兔抗 YPEL1 多克隆抗体

中文名称：兔抗 YPEL1 多克隆抗体

英文名称：Anti-YPEL1 rabbit polyclonal antibody

别名：FKSG3

抗原：YPEL1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is located in the region associated with DiGeorge syndrome on chromosome 22. The encoded protein localizes to the centrosome and nucleolus and may play a role in the regulation of cell division.
Applications:	ELISA, IHC
Name of antibody:	YPEL1
Immunogen:	Full length fusion protein
Full name:	yippee like 1
Synonyms:	FKSG3
SwissProt:	O60688
ELISA Recommended diluti	5000-10000

on:	
IHC positive control:	Human colorectal cancer and human cervical cancer
IHC Recommend dilution:	25-100



