

兔抗 YPEL2 多克隆抗体

中文名称：兔抗 YPEL2 多克隆抗体

英文名称：Anti-YPEL2 rabbit polyclonal antibody

别名：FKSG4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：YPEL2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

YPEL2 (yippee-like 2) belongs to a family of five yippee-like proteins, all of which localize to the centrosome or mitotic spindle and are widely expressed in both adult and fetal tissue. This localization plus the fact that the family of human YPEL proteins share a high degree of sequence homology across species suggests that these proteins may have a conserved function involved in cell division. YPEL2 might be an important factor during the development and malignant transformation of tissues, most notably pancreatic and breast tumors.

Applications:	ELISA, IHC
Name of antibody:	YPEL2
Immunogen:	Full length fusion protein
Full name:	yippee like 2
Synonyms:	FKSG4
SwissProt:	Q96QA6
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and human thyroid cancer
IHC Recommend dilution:	25-100



