

兔抗 YPEL4 多克隆抗体

中文名称：兔抗 YPEL4 多克隆抗体

英文名称：Anti-YPEL4 rabbit polyclonal antibody

抗原：YPEL4

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	YPEL4 (yippee-like 4) belongs to a family of five yippee-like proteins, all of which localize to the centrosome or mitotic spindle and are widely expressed in both adult and fetal tissue. This localization plus the fact that the family of human YPEL proteins share a high degree of sequence homology across species suggests that these proteins may have a conserved function involved in cell division. YPEL4 is ubiquitously expressed in adult tissues and has been shown to associate with the major vault protein (MVP). It has been suggested that MVP can inhibit YPEL4's ability to activate Elk-1 in the MAPK signaling pathway.
Applications:	ELISA, IHC

Name of antibody:	YPEL4
Immunogen:	Synthetic peptide of human YPEL4
Full name:	yippee like 4
SwissProt:	Q96NS1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and human lung cancer
IHC Recommend dilution:	25-100



