

兔抗 KDF1 多克隆抗体

中文名称: 兔抗 KDF1 多克隆抗体

英文名称: Anti-KDF1 rabbit polyclonal antibody

别名: C1orf172

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: KDF1

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Plays a role in the regulation of the epidermis formation during early development. Required both as an inhibitor of basal cell proliferation and a promoter of differentiation of basal progenitor cell progeny (By similarity).
Applications:	ELISA, IHC
Name of antibody:	KDF1
Immunogen:	Synthetic peptide of human KDF1
Full name:	keratinocyte differentiation factor 1
Synonyms:	C1orf172

SwissProt:	Q8NAX2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

