

兔抗 YPEL5 多克隆抗体

中文名称：兔抗 YPEL5 多克隆抗体

英文名称：Anti-YPEL5 rabbit polyclonal antibody

级别：CGI-127

抗原：YPEL5

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	YPEL5 (yippee-like 5), also known as CGI-127, is a 121 amino acid protein that may be involved in cell division-related function. During cell cycle progression, YPEL5 localizes to multiple subcellular regions. At interphase of mitosis, YPEL5 localizes to the nucleus and centrosome, then changes its location sequentially to the spindle poles, mitotic spindle and spindle midzone, and finally it is transferred to the midbody at cytokinesis. The function of YPEL5 during cell division is not yet fully understood.
Applications:	ELISA, IHC

Name of antibody:	YPEL5
Immunogen:	Synthetic peptide of human YPEL5
Full name:	yippee like 5
Synonyms:	CGI-127
SwissProt:	P62699
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain and human esophagus cancer
IHC Recommend dilution:	40-250



