

兔抗 YWHAB 多克隆抗体

中文名称：兔抗 YWHAB 多克隆抗体

英文名称： Anti-YWHAB rabbit polyclonal antibody

别名： HS1, GW128, YWHAA, KCIP-1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： YWHAB

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein belonging to the 14-3-3 family of proteins, members of which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals. The encoded protein has been shown to interact with RAF1 and CDC25 phosphatases, suggesting that it may play a role in linking mitogenic signaling and the cell cycle machinery. Two transcript variants, which encode the same protein, have been identified for this gene.
Applications:	ELISA, WB
Name of antibody:	YWHAB

Immunogen:	Synthetic peptide of human YWHAB
Full name:	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, beta polypeptide
Synonyms :	HS1, GW128, YWHAA, KCIP-1
SwissProt:	P31946
ELISA Recommended dilution:	1000-5000
WB Predicted band size:	28 kDa
WB Positive control:	Hela cells
WB Recommended dilution:	500-2000

