

兔抗 SLC2A6 多克隆抗体

中文名称：兔抗 SLC2A6 多克隆抗体

英文名称： Anti-SLC2A6 rabbit polyclonal antibody

别名： GLUT6; GLUT9; HSA011372

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC2A6

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Members of the Glut family may enhance the metabolic activity of tumor cells. Glut6 is part of the third out of three classes of Gluts. Glut6 is mainly expressed in the brain, spleen and peripheral leukocytes. It appears to be regulated by subcellular redistribution, because it is targeted to intracellular compartments by di-leucine motifs, recycling itself in a Dynamin-dependent manner.
Applications:	ELISA, WB, IHC
Name of antibody:	SLC2A6
Immunogen:	Synthetic peptide of human SLC2A6
Full name:	solute carrier family 2 (facilitated glucose transporter), me

	mber 6
Synonyms :	GLUT6; GLUT9; HSA011372
SwissProt:	Q9UGQ3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human prostate cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	54 kDa
WB Positive control:	Human fetal liver tissue
WB Recommended dilution:	200-1000



