

兔抗 KDM2B 多克隆抗体

中文名称：兔抗 KDM2B 多克隆抗体

英文名称：Anti-KDM2B rabbit polyclonal antibody

别名：lysine demethylase 2B; CXXC2; Fbl10; PCCX2; FBXL10; JHDM1B

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：KDM2B

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbls class. Multiple alternatively spliced tra

	nscript variants have been found for this gene, but the full-length nature of some variants has not been determined.
Applications:	ELISA, IHC
Name of antibody:	KDM2B
Immunogen:	Synthetic peptide of human KDM2B
Full name:	lysine demethylase 2B
Synonyms:	CXXC2; Fbl10; PCCX2; FBXL10; JHDM1B
SwissProt:	Q8NHM5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer and Human brain
IHC Recommend dilution:	50-200



