

兔抗 YWHAH 多克隆抗体

中文名称：兔抗 YWHAH 多克隆抗体

英文名称： Anti-YWHAH rabbit polyclonal antibody

别名： YWHA1

抗原： YWHAH

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and bovine orthologs. This gene contains a 7 bp repeat sequence in its 5' UTR, and changes in the number of this repeat have been associated with early-onset schizophrenia and psychotic bipolar disorder
Applications:	WB, IF

Name of antibody:	YWHAH
Immunogen:	Synthesized peptide derived from internal of human 14-3-3 η.
Full name:	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, eta
Synonyms :	YWHA1
SwissProt:	Q04917
WB Predicted band size:	28 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000
IF positive control:	Hela cells
IF Recommend dilution:	100-500



