

兔抗 PLCL2 多克隆抗体

中文名称：兔抗 PLCL2 多克隆抗体

英文名称：Anti-PLCL2 rabbit polyclonal antibody

别名：phospholipase C like 2; PLCE2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PLCL2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May play an role in the regulation of Ins(1,4,5)P3 around the endoplasmic reticulum.
Applications:	ELISA, IHC
Name of antibody:	PLCL2
Immunogen:	Fusion protein of human PLCL2
Full name:	phospholipase C like 2
Synonyms:	PLCE2
SwissProt:	Q9UPR0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human tonsil
IHC Recommend dilution:	25-100

