

兔抗 KDM3B 多克隆抗体

- 中文名称：兔抗 KDM3B 多克隆抗体
- 英文名称：Anti-KDM3B rabbit polyclonal antibody
- 别名：lysine demethylase 3B; 5qNCA; NET22; C5orf7; JMJD1B
- 相关类别：一抗
- 抗原：KDM3B
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	Histone demethylase that specifically demethylates 'Lys-9' of histone H3, thereby playing a central role in histone code. Demethylation of Lys residue generates formaldehyde and succinate. May have tumor suppressor activity.
Applications:	ELISA, IHC
Name of antibody:	KDM3B
Immunogen:	Synthetic peptide of human KDM3B
Full name:	lysine demethylase 3B
Synonyms:	5qNCA; NET22; C5orf7; JMJD1B

SwissProt:	Q7LBC6
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	30-150

