

## 兔抗 SLC2A12 多克隆抗体

中文名称：兔抗 SLC2A12 多克隆抗体

英文名称： Anti-SLC2A12 rabbit polyclonal antibody

别名： solute carrier family 2 member 12; GLUT8; GLUT12

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC2A12

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	SLC2A12 belongs to a family of transporters that catalyze the uptake of sugars through facilitated diffusion (Rogers et al., 2002). This family of transporters show conservation of 12 transmembrane helices as well as functionally significant amino acid residues (Joost and Thorens, 2001 [PubMed 11780753]).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SLC2A12
<b>Immunogen:</b>	Synthetic peptide of human SLC2A12
<b>Full name:</b>	solute carrier family 2 member 12
<b>Synonyms:</b>	GLUT8; GLUT12

<b>SwissProt:</b>	Q8TD20
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	50-100

