

兔抗 SLC4A1AP 多克隆抗体

中文名称：兔抗 SLC4A1AP 多克隆抗体

英文名称： Anti-SLC4A1AP rabbit polyclonal antibody

别名： HLC3

抗原： SLC4A1AP

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Solute carrier family 4 anion exchanger member 1 adapter protein (SLC4A1AP), also designated Kanadaptin or human lung cancer oncogene 3 protein (HLC-3), is a 796 amino acid protein that is widely expressed in many tissues, including kidney, lung, liver, brain and skeletal and cardiac muscle. SLC4A1AP is a multidomain protein that localizes to the nucleus where it may play a role in signaling. SLC4A1AP was previously thought to act as an adaptor protein or chaperone involved in targeting kAE1 to the plasma membrane. However, recent studies suggest SLC4A1AP does not interact with kA

	E1. The gene encoding SLC4A1AP maps to chromosome 2, which consists of 237 million bases and makes up approximately 8% of the human genome. A number of genetic diseases are linked to genes on chromosome 2, including Harlequin ichthyosis, sitosterolemia and Alström syndrome.
Applications:	WB
Name of antibody:	SLC4A1AP
Immunogen:	Synthesized peptide derived from internal of human NADAP.
Full name:	solute carrier family 4 (anion exchanger), member 1, adaptor protein
Synonyms:	HLC3
SwissProt:	Q9BWU0
WB Predicted band size:	89 kDa
WB Positive control:	Huvec cells lysate
WB Recommended dilution:	500-3000

