

## 兔抗 YWHAZ (Phospho-Ser58)多克隆抗体

中文名称: 兔抗 YWHAZ (Phospho-Ser58)多克隆抗体

英文名称: Anti-YWHAZ (Phospho-Ser58) rabbit polyclonal antibody

别 名: YWHAD; KCIP-1; 14-3-3-zeta

相关类别: 一抗

储 存: 冷冻(-20℃) 避光

宿 主: Rabbit

抗 原: YWHAZ (Phospho-Ser58)

反应种属: Human, Mouse, Rat

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

## 技术规格

Background:	This gene product belongs to the 14-3-3 family of prote ins which mediate signal transduction by binding to pho sphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and shee porthologs. The encoded protein interacts with IRS1 protein, suggesting a role in regulating insulin sensitivity. Several transcript variants that differ in the 5' UTR but that encode the same protein have been identified for this gene.
Applications:	WB, IHC, IF



Name of antibody:	YWHAZ (Phospho-Ser58)
Immunogen:	Synthetic peptide of human YWHAZ (Phospho-Ser58)
Full name:	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, zeta polypeptide (Phospho-Ser58)
Synonyms:	YWHAD; KCIP-1; 14-3-3-zeta
SwissProt:	P63104
IHC positive control:	Human breast carcinoma
IHC Recommend dilution:	50-100
WB Predicted band size:	28 kDa
WB Positive control:	Hela cells treated with UV
WB Recommended dilution:	500-1000
IF Positive control:	Hela cells
IF Recommended dilution	100-200









