

兔抗 KDM5A 多克隆抗体

中文名称: 兔抗 KDM5A 多克隆抗体

英文名称: Anti-KDM5A rabbit polyclonal antibody

别名: RBP2; RBBP2; RBBP-2

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: KDM5A

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the Jumonji, AT-rich interactive domain 1 (JARID1) histone demethylase protein family . The encoded protein plays a role in gene regulation through the histone code by specifically demethylating lysine 4 of histone H3. The encoded protein interacts with many other proteins, including retinoblastoma protein, and is implicated in the transcriptional regulation of Hox genes and cytokines. This gene may play a role in tumor progression.
Applications:	ELISA, WB, IHC
Name of antibody:	KDM5A

Immunogen:	Synthetic peptide of human KDM5A
Full name:	lysine (K)-specific demethylase 5A
Synonyms :	RBP2; RBBP2; RBBP-2
SwissProt:	P29375
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	192 kDa
WB Positive control:	231 and hela cells
WB Recommended dilution:	500-2000



