

兔抗 PLIN3 多克隆抗体

中文名称：兔抗 PLIN3 多克隆抗体

英文名称： Anti-PLIN3 rabbit polyclonal antibody

别名： PP17; TIP47; M6PRBP1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PLIN3

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Mannose 6-phosphate receptors (MPRs) deliver lysosomal hydrolase from the Golgi to endosomes and then return to the Golgi complex. The protein encoded by this gene interacts with the cytoplasmic domains of both cation-independent and cation-dependent MPRs, and is required for endosome-to-Golgi transport. This protein also binds directly to the GTPase RAB9 (RAB9A), a member of the RAS oncogene family. The interaction with RAB9 has been shown to increase the affinity of this protein for its cargo. Multiple transcript variants encoding different isoforms have been found for this gene.

Applications:	ELISA, IHC
Name of antibody:	PLIN3
Immunogen:	Synthetic peptide of human PLIN3
Full name:	perilipin 3
Synonyms :	PP17; TIP47; M6PRBP1
SwissProt:	O60664
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human brain
IHC Recommend dilution:	50-200

