

## 兔抗 PLIN4 多克隆抗体

中文名称：兔抗 PLIN4 多克隆抗体

英文名称： Anti-PLIN4 rabbit polyclonal antibody

别名： perilipin 4; S3-12; KIAA1881

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PLIN4

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Members of the perilipin family, such as PLIN4, coat intracellular lipid storage droplets (Wolins et al., 2003 [PubMed 12840023]). May play a role in triacylglycerol packaging into adipocytes. May function as a coat protein involved in the biogenesis of lipid droplets (By similarity).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PLIN4
<b>Immunogen:</b>	Synthetic peptide of human PLIN4
<b>Full name:</b>	perilipin 4
<b>Synonyms:</b>	S3-12; KIAA1881
<b>SwissProt:</b>	Q96Q06

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human prostate cancer
<b>IHC Recommend dilution:</b>	50-200



