

## 兔抗 SLC5A11 多克隆抗体

中文名称：兔抗 SLC5A11 多克隆抗体

英文名称： Anti-SLC5A11 rabbit polyclonal antibody

别名： KST1; RKST1; SGLT6; SMIT2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC5A11

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Cotransporters, such as SLC5A11, represent a major class of proteins that make use of ion gradients to drive active transport for the cellular accumulation of nutrients, neurotransmitters, osmolytes, and ions Roll et al.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SLC5A11
<b>Immunogen:</b>	Synthetic peptide of human SLC5A11
<b>Full name:</b>	solute carrier family 5 (sodium/inositol cotransporter), member 11
<b>Synonyms :</b>	KST1; RKST1; SGLT6; SMIT2

<b>SwissProt:</b>	Q8WWX8
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	30-150

