

兔抗 KHSRP 多克隆抗体

中文名称:	兔抗 KHSRP 多克隆抗体
英文名称:	Anti-KHSRP rabbit polyclonal antibody
别名:	KH-type splicing regulatory protein; FBP2; KSRP; FUBP2
相关类别:	一抗
储存:	冷冻 (-20℃)
宿主:	Rabbit
抗原:	KHSRP
反应种属:	Human, Mouse, Rat
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	The KHSRP gene encodes a multifunctional RNA-binding protein implicated in a variety of cellular processes, including transcription, alternative pre-mRNA splicing, and mRNA localization.
Applications:	ELISA, WB, IHC
Name of antibody:	KHSRP
Immunogen:	Synthetic peptide of human KHSRP
Full name:	KH-type splicing regulatory protein
Synonyms:	FBP2; KSRP; FUBP2
SwissProt:	Q92945
ELISA Recommended dilutio	5000-10000

n:	
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	73 kDa
WB Positive control:	HepG2 and Hela cell lysates
WB Recommended dilution:	1000-5000



