

兔抗 ZAR1 多克隆抗体

中文名称：兔抗 ZAR1 多克隆抗体

英文名称：Anti-ZAR1 rabbit polyclonal antibody

别名：Z3CXXC6

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ZAR1

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This maternal effect gene is oocyte-specific and encodes a protein that is thought to function in the initiation of embryogenesis. A similar protein in mouse is required for female fertility.
Applications:	WB
Name of antibody:	ZAR1
Immunogen:	Synthesized peptide derived from C-terminal of human ZAR1.
Full name:	zygote arrest 1
Synonyms :	Z3CXXC6
SwissProt:	Q86SH2

WB Predicted band size:	46 kDa
WB Positive control:	293 cells lysate
WB Recommended dilution:	500-3000

